



# Center for Professional Practice of Nursing

COURSE NAME: Trainer Training the Trainer POCT

COURSE NUMBER: DAHS-NGNPCTT23-POCT

DEPARTMENT NAME: \_\_\_\_\_

2023 POCT	NAME	Net ID	Date	Checklist Nova StatStrip	Checklist Coloscreen	Checklist Gastrocult	Checklist Urine Dipstick	Checklist Urine Pregnancy	Checklist Bacterial Vaginosis	Checklist CoaguCheck	Checklist Cobas Liat	Checklist HemoCue Hb 201 DM	Checklist pH paper	Checklist Alere iCup	Hemoglobin A1c (HbA1c) Checklist	Adeno Plus	Initials	
	1				<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	
	2				<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	
	3				<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	
	4				<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	
	5				<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	
	6				<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	
	7				<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	
	8				<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	
	9				<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	
10				<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>		
11				<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>		

## 2023 POINT OF CARE TESTING

<b>Name:</b>	<b>Employee ID#:</b>
<b>Unit:</b>	<b>Title:</b>
<b>Due Date:</b>	
<b>Scan Page 1 &amp; 2 and email to: <a href="mailto:hs-cppn@ucdavis.edu">hs-cppn@ucdavis.edu</a></b> Preceptor or other verified personnel date and initial completed skill.	

SIGNATURE PAGE:		
Signature and Printed Name of Verifier (preceptor or other verified personnel) who have initialed on this form:		
Initial:	Print Name:	Signature:

**PRECEPTEE STATEMENT AND SIGNATURE:**

I have read and understand the appropriate UCDH Policies and Procedures and/or equipment operations manual, I have demonstrated the ability to perform the verified skills asnoted, and I have the knowledge of the resources available to answer questions.

<b>Printed Name</b>	<b>Signature</b>

**Scan Page 1 and 2 and email to: [hs-cppn@ucdavis.edu](mailto:hs-cppn@ucdavis.edu)**

**2023 POINT OF CARE TESTING**

<b>Name:</b>	<b>Employee ID#:</b>
<b>Unit:</b>	<b>Title:</b>

**These skills will be considered complete when all below performance criteria are completed.**

AdenoPlus by RPS # DAHS-NSCAPRPS23-POCT	Date	Verifier Initials
<p><b>REFERENCES:</b></p> <ol style="list-style-type: none"> <li><a href="#">UC Davis Health Policy POCT.57: AdenoPlus by RPS</a></li> <li>RPS Package Insert AdenoPlus External Controls. Rapid Pathogen Screening Inc., Sarasota ,FL 34240</li> </ol>		
<ul style="list-style-type: none"> <li><input type="checkbox"/> Completes online module and reads associated policy</li> <li><input type="checkbox"/> Describes sample collection and handling.</li> <li><input type="checkbox"/> Demonstrates proper QC and patient testing procedure of AdenoPlus test</li> <li><input type="checkbox"/> Correctly interprets patient and quality control results</li> <li><input type="checkbox"/> Verifies test kit is functioning properly, supplies and reagents are not outdated, and quality controls are within acceptable limits</li> <li><input type="checkbox"/> Describe Quality Control Remedial Action for out-of-control results</li> </ul>		

**DO NOT SEND CHECKLIST TO CPPN**

**File Checklists in Home Department for 2 Years**

**2023 POINT OF CARE TESTING**

Name:	Employee ID#:
Unit:	Title:
<b>These skills will be considered complete when all below performance criteria are completed.</b>	

<b>Bacterial Vaginosis Rapid Test Checklist 2023</b>	<b># DAHS-NSCBVRT23-POCT</b>		
<b>REFERENCES:</b>		<b>Date</b>	<b>Verifier Initials</b>
1. <a href="#">UC Davis Health Policy POCT.58, Bacterial Vaginosis Rapid Test by OSOM BV Blue Test, Sekisui Diagnostics</a>			
<ul style="list-style-type: none"> <li><input type="checkbox"/> Completes online module and reads associated policy</li> <li><input type="checkbox"/> Describes proper sample collection of vaginal fluid</li> <li><input type="checkbox"/> Demonstrates proper test procedure</li> <li><input type="checkbox"/> Correctly interprets results</li> <li><input type="checkbox"/> Verbalizes understanding of valid results (internal procedural controls). Articulates appropriate action for out-of-control results</li> <li><input type="checkbox"/> Demonstrates quality control testing using Sekisui OSOM BVBLUE positive and negative sialidase controls. Describes stability of controls and when testing must be performed</li> </ul>			

**DO NOT SEND CHECKLIST TO CPPN**

**File Checklists in Home Department for 2 Years**

## 2023 POINT OF CARE TESTING

Name:	Employee ID#:
Unit:	Title:
<b>These skills will be considered complete when all below performance criteria are completed.</b>	

CoaguChek Checklist 2023	# DAHS-NSCCOAGU23-POCT		
REFERENCES:		Date	Verifier Initials
1. <a href="#">UC Davis Health Policy POCT.42: Prothrombin Time (INR), Whole Blood by Roche CoaguChek XS</a>			
<ul style="list-style-type: none"> <li><input type="checkbox"/> Completes online module and reads associated policy</li> <li><input type="checkbox"/> Describe the Roche CoaguChek XS Meter components and their function</li> <li><input type="checkbox"/> Describe the indications for use</li> <li><input type="checkbox"/> Describe specimen collection criteria and limitations</li> <li><input type="checkbox"/> Describe storage requirements, stability criteria, and handling of testing strips</li> <li><input type="checkbox"/> Demonstrate how to ensure matching of meter and test strips with the code chips</li> <li><input type="checkbox"/> Describe the integrated quality controls</li> <li><input type="checkbox"/> Describe remedial action for Out-of-Control Results</li> </ul>			

**DO NOT SEND CHECKLIST TO CPPN**

**File Checklists in Home Department for 2 Years**

**2023 POINT OF CARE TESTING**

<b>Name:</b>	<b>Employee ID#:</b>
<b>Unit:</b>	<b>Title:</b>

**These skills will be considered complete when all below performance criteria are completed.**

<b>Cobas® Liat® System Checklist 2023 # DAHS-NSCCOBAS23-POCT</b>		<b>Date</b>	<b>Verifier Initials</b>
<b>REFERENCES:</b> 1. UC Davis Health Dept. Pathology & Laboratory Medicine (QC book located in the ED) 2. <a href="#">UC Davis Health Policy POCT.59: Roche cobas Liat System, Molecular testing of Infectious Disease</a>			
<input type="checkbox"/> Completes online module and reads associated policy <input type="checkbox"/> Describes sample collection and handling <input type="checkbox"/> Describe the stability and handling of the reagent kits <input type="checkbox"/> Verbalizes understanding of valid results. Articulates what to do when an indeterminate or invalid result occurs <input type="checkbox"/> Direct observations of routine patient testing, including, as applicable, patient identification and sample processing and testing <input type="checkbox"/> Demonstrates routine cleaning of instrumentation and immediate work area			

**DO NOT SEND CHECKLIST TO CPPN**

**File Checklists in Home Department for 2 Years**

**2023 POINT OF CARE TESTING**

<b>Name:</b>	<b>Employee ID#:</b>
<b>Unit:</b>	<b>Title:</b>
<b>These skills will be considered complete when all below performance criteria are completed.</b>	

<b>ColoScreen Checklist 2023</b>	<b># DAHS-NSCCOLO23-POCT</b>		
<b>REFERENCES:</b>		<b>Date</b>	<b>Verifier Initials</b>
1. <a href="#">UC Davis Health Policy POCT.19: Occult Blood, Fecal by Helena Laboratories ColoScreen</a> 2. ColoScreen: A Test for Occult Blood, Package Insert, Helena Laboratories			
<ul style="list-style-type: none"> <li><input type="checkbox"/> Completes online module and associated policy</li> <li><input type="checkbox"/> Identify factors that may affect the guaiac test</li> <li><input type="checkbox"/> Describes collection of the stool specimen (able to teach patient for home use)</li> <li><input type="checkbox"/> Demonstrates proper testing procedure using ColoScreen developer</li> <li><input type="checkbox"/> Correctly interprets specimen results</li> <li><input type="checkbox"/> Demonstrates how to perform QC using the internal ColoCheck Monitors. Verbalizes understanding of procedure to follow when control monitors fail</li> <li><input type="checkbox"/> Describes the stability and handling of the guaiac slides and ColoScreen developer</li> </ul>			

**DO NOT SEND CHECKLIST TO CPPN**

**File Checklists in Home Department for 2 Years**

**2023 POINT OF CARE TESTING**

<b>Name:</b>	<b>Employee ID#:</b>
<b>Unit:</b>	<b>Title:</b>

**These skills will be considered complete when all below performance criteria are completed.**

<b>Gastrocult Checklist 2023</b>	<b># DAHS-NSCGAST23-POCT</b>		
<b>REFERENCES:</b> 1. <a href="#">UC Davis Health Policy POCT.20: Occult Blood, Gastric Fluid</a> 2. Gastrocult, Package Insert, Beckman Coulter		<b>Date</b>	<b>Verifier Initials</b>
<ul style="list-style-type: none"> <li><input type="checkbox"/> Completes online module and reads associated policy</li> <li><input type="checkbox"/> Identify factors that may affect the Gastrocult test</li> <li><input type="checkbox"/> Demonstrates proper testing procedure of gastric aspirate, using Gastrocult developer</li> <li><input type="checkbox"/> Correctly interprets results</li> <li><input type="checkbox"/> Demonstrates how to perform QC using internal Performance Monitors. Verbalizes understanding of procedure to follow when control monitors fail</li> <li><input type="checkbox"/> Describes the stability and handling of the Gastrocult slides and developer</li> </ul>			

**DO NOT SEND CHECKLIST TO CPPN**

**File Checklists in Home Department for 2 Years**



## 2023 POINT OF CARE TESTING

<b>Name:</b>	<b>Employee ID#:</b>
<b>Unit:</b>	<b>Title:</b>

**These skills will be considered complete when all below performance criteria are completed.**

HemoCue Hb 201 DM Checklist 2023	# DAHS-NSCHEMOC23-POCT		
<b>REFERENCES:</b> 1. <a href="#">UC Davis Health Policy POCT.43: Hemoglobin, Whole Blood by HemoCue Hb 210 DM</a> 2. HemoCue Hb 201 DM Analyzer Instructions for Use Manual, HemoCue 3. HemoCue Hb 201 DM Analyzer & HemoCue DM Docking Station Reference Manual, HemoCue 4. HemoCue Hb 201 Microcuvettes Package Insert, HemoCue 5. HemoCue Cleaner Package Insert, HemoCue 6. HemoCue Technical Letter No 2 September 2010, HemoCueAB 7. R & D Systems R & D Glu/Hgb Control Package Insert, 6/07, R & D Systems, Inc.		<b>Date</b>	<b>Verifier Initials</b>
<input type="checkbox"/> Completes online module and reads associated policy <input type="checkbox"/> Describe how to handle the microcuvettes and proper sampling technique <input type="checkbox"/> Describe the storage and stability of the liquid controls <input type="checkbox"/> Demonstrate how to perform daily QC using liquid controls <input type="checkbox"/> Describe Quality Control Remedial Action for out-of-control results <input type="checkbox"/> Describe how to perform patient testing on the HemoCue Hb 201 DM Analyzer <input type="checkbox"/> Describe what to do when you have a critical value <input type="checkbox"/> Describe the limitations of the HemoCue HB 201 DM procedure <input type="checkbox"/> Describe the maintenance for the HemoCue Hb 201 DM Analyzer			

**DO NOT SEND CHECKLIST TO CPPN**

**File Checklists in Home Department for 2 Years**

## 2023 POINT OF CARE TESTING

<b>Name:</b>	<b>Employee ID#:</b>
<b>Unit:</b>	<b>Title:</b>

**These skills will be considered complete when all below performance criteria are completed.**

Hemoglobin A1c Checklist 2023	# DAHS-NSCHEMOG23-POCT	
<b>REFERENCES:</b> 1. <a href="#">UC Davis Health Policy POCT.40: Hemoglobin A1c, Whole Blood by Siemens DCA Vantage Analyzer</a> 2. Siemens DCA Systems Hemoglobin A1c Normal and Abnormal Control Kit Package Insert 3. Siemens DCA Systems Hemoglobin A1c Reagent Kit Package Insert 4. Siemens DCA Vantage Analyzer Operator's Guide 5. Siemens DCA Vantage Quick Set-up Guide		<b>Verifier Initials</b>
<ul style="list-style-type: none"> <li><input type="checkbox"/> Completes online module and reads associated policy</li> <li><input type="checkbox"/> Verifies physician written order or standardized procedure</li> <li><input type="checkbox"/> Identifies the patient using two forms of identifiers</li> <li><input type="checkbox"/> Explains the procedure and reason for the test to the patient/family</li> <li><input type="checkbox"/> Gathers supplies/equipment</li> <li><input type="checkbox"/> Verifies test kit/machine/equipment is functioning properly, supplies and reagents are not outdated, and quality controls are within acceptable limits.</li> <li><input type="checkbox"/> Appropriately documents result on correct form</li> <li><input type="checkbox"/> Notifies provider of abnormal results</li> <li><input type="checkbox"/> Appropriately charges for supplies and procedure when applicable</li> </ul>	<b>Date</b>	

**DO NOT SEND CHECKLIST TO CPPN**

**File Checklists in Home Department for 2 Years**

## 2023 POINT OF CARE TESTING

<b>Name:</b>	<b>Employee ID#:</b>
<b>Unit:</b>	<b>Title:</b>
<b>These skills will be considered complete when all below performance criteria are completed</b>	

Nova StatStrip Checklist 2023	# DAHS-NSCNOVA23-POCT		
References:		Date	Verifier Initials
<ol style="list-style-type: none"> <li>1. Instructions for Use: Nova Biomedical StatStrip Glucose Hospital Meter System. Version 1.86. UC Davis Health</li> <li>2. <a href="#">UC Davis Health Policy POCT.51: Glucose, Whole Blood by Nova StatStrip Glucose Hospital Meter System</a></li> </ol>			
<ul style="list-style-type: none"> <li><input type="checkbox"/> Completes online module and reads associated policy</li> <li><input type="checkbox"/> Describes the Nova StatStrip Glucose Hospital Meter System components and their functions</li> <li><input type="checkbox"/> Describes the stability and handling of glucose test strips and controls, including open dating, expiration dates and procedure notes</li> <li><input type="checkbox"/> Describes specimen requirements: amount, types, and recommended anticoagulants</li> <li><input type="checkbox"/> Describes the conditions when it is not appropriate to perform a glucose fingerstick</li> <li><input type="checkbox"/> Demonstrates how to perform quality control testing using one control solution. Describes when QC must be run, what QC Lockout is, and what to do if the QC test results FAIL</li> <li><input type="checkbox"/> Demonstrate how to run a patient test. (Using one control run as a patient)</li> <li><input type="checkbox"/> Describes the procedure for critical values and questionable patient test results</li> <li><input type="checkbox"/> Describes the limitations of the Nova StatStrip Glucose Hospital Meter System glucose test results</li> <li><input type="checkbox"/> Demonstrates how to review results on the meter</li> <li><input type="checkbox"/> Describes the routine maintenance for the Nova StatStrip Glucose Hospital Meter System, docking station, and accessory box</li> </ul>			

**DO NOT SEND CHECKLIST TO CPPN**

**File Checklists in Home Department for 2 Years**

**2023 POINT OF CARE TESTING**

Name:	Employee ID#:
Unit:	Title:

**These skills will be considered complete when all below performance criteria are completed.**

pH Paper Checklist 2023	# DAHS-NSCPHP23-POCT	
<b>REFERENCES:</b> 1. <a href="#">UC Davis Health Policy POCT.23: pH Fluid Nitrazine Paper</a> 2. <a href="#">UC Davis Health Policy POCY.23: pH Fluid by pHDracid Paper</a> 3. Instructions for Use: pHizatest Product Insert. Rev. Nov 15, 2019 4. pHizatest Specification Sheet, version #4, 06-19-2013. 5. Micro Essentials pHizatest Product Information Sheet, Rev. 12-20-12. 6. Micro Essential Hydrion Product Information Sheet, Rev. 10-16-13.	Date	Verifier Initials
<ul style="list-style-type: none"> <li><input type="checkbox"/> Completes online module and reads associated policies</li> <li><input type="checkbox"/> Describes indications of use for each paper type</li> <li><input type="checkbox"/> Describes the stability and handling, including open dating, expiration dates and quality control</li> <li><input type="checkbox"/> Describes/demonstrates how to perform and result a test</li> <li><input type="checkbox"/> Describes the limitations of the procedure</li> </ul>		

**DO NOT SEND CHECKLIST TO CPPN**

**File Checklists in Home Department for 2 Years**

**2023 POINT OF CARE TESTING**

<b>Name:</b>	<b>Employee ID#:</b>
<b>Unit:</b>	<b>Title:</b>

**These skills will be considered complete when all below performance criteria are completed.**

Urine Dipstick/Uristix Checklist 2023	# DAHS-NSCURID23-POCT	
<p><b>REFERENCES:</b></p> <ol style="list-style-type: none"> <li><a href="#">UC Davis Health Policy: POCT.54: Urine Dipstick by Siemens Uristix and Multistix 10SG Reagent Strips</a></li> <li>MAS UA Control, Package insert, Thermo Scientific.</li> <li>Multistix 10 SG Reagent Strips, Package insert, Siemens.</li> <li>Uristix Reagent Strips, package insert, Siemens.</li> <li>Keto-Diastix Reagent Strips, Siemens.</li> </ol>		<b>Verifier Initials</b>
<ul style="list-style-type: none"> <li><input type="checkbox"/> Completes online module and reads associated policy</li> <li><input type="checkbox"/> Collects appropriate urine specimen in clean, dry container</li> <li><input type="checkbox"/> Demonstrates proper test procedure</li> <li><input type="checkbox"/> Properly compares the reagent areas to the corresponding color chart on the bottle label at the times specified</li> <li><input type="checkbox"/> Demonstrates quality control testing using the MAS UA controls. Describes stability of controls and when controls must be run</li> <li><input type="checkbox"/> Verbalizes understanding of procedure to follow when control results are out-of-range</li> </ul>	<b>Date</b>	

**DO NOT SEND CHECKLIST TO CPPN**

**File Checklists in Home Department for 2 Years**

**2023 POINT OF CARE TESTING**

Name:	Employee ID#:
Unit:	Title:
<b>These skills will be considered complete when all below performance criteria are completed</b>	

Urine Drug Screen by Alere iCup® 2023	# DAHS-NSCALERE23-POCT	Date	Verifier Initials
<b>REFERENCES:</b> 1. <a href="#">UC Davis Health Policy POCT.60: Urine Drug Screen by Alere iCup DX Pro 2</a>			
<input type="checkbox"/> Completes online module and reads associated policy <input type="checkbox"/> Describes proper sample collection of urine and demonstrates importance of appropriate collection times <input type="checkbox"/> Demonstrates or describes proper test procedure <input type="checkbox"/> Correctly interprets results <input type="checkbox"/> Verbalizes understanding of valid results (internal procedural controls) and articulates appropriate action for out-of-control results <input type="checkbox"/> Describes external quality control testing procedure, stability of reagents and requirements of performance			

**DO NOT SEND CHECKLIST TO CPPN**

**File Checklists in Home Department for 2 Years**

**2023 POINT OF CARE TESTING**

<b>Name:</b>	<b>Employee ID#:</b>
<b>Unit:</b>	<b>Title:</b>

**These skills will be considered complete when all below performance criteria are completed.**

Urine Pregnancy Test Checklist 2023	# DAHS-NSCURIP23-POCT	
<b>REFERENCES:</b> 1. <a href="#">UC Davis Health Policy POCT.41: hCG Pregnancy Test, Urine by Sure-Vue hCG</a> 2. Sure-Vue® Urine hCG, Package Insert, Fisher Scientific Company 3. Sure-Vue® hCG Urine Control Set, Package Insert, Fisher Scientific Company	Date	Verifier Initials
<ul style="list-style-type: none"> <li><input type="checkbox"/> Completes online module and reads associated policy</li> <li><input type="checkbox"/> Collects appropriate urine specimen in clean, dry container. Equilibrate urine and test device to room temperature</li> <li><input type="checkbox"/> Demonstrates proper test procedure</li> <li><input type="checkbox"/> Correctly interprets results</li> <li><input type="checkbox"/> Verbalizes understanding of valid results (internal procedural controls). Articulates what to do when there is no control line visible</li> <li><input type="checkbox"/> Demonstrates quality control testing using Sure-Vue positive and negative hCG controls. Describes stability of controls and when they must be run</li> </ul>		

**DO NOT SEND CHECKLIST TO CPPN**

**File Checklists in Home Department for 2 Years**