

# Comfort Commitment: ShotBlocker® Use Guide

The ShotBlocker® uses blunt contact points to distract patients 2 years and older from the potential pain of needle procedures.

## Overview:

- Select injection site and prep skin per protocol.
- Hold the ShotBlocker® so that the blunt contact points touch the patient’s skin at the injection site.
- Press ShotBlocker® firmly against the skin. Do not move the ShotBlocker® until the injection is complete.
- Administer injection through or near the ShotBlocker® opening. If more than 20 seconds lapse between the placement of ShotBlocker® and the procedure, remove the ShotBlocker® from the skin and repeat the process.
- The ShotBlocker® is for single-patient use. Patients can keep for future use. Wash with soap and water between uses.

## When introducing to patients:

- Explain how the blunt contact points work to block the pain from the injection.
- Let the patient feel the tool and show how it will be used.



**Product information:** [Bionix.com/shotblocker.html](https://bionix.com/shotblocker.html)

**Website:** <https://ucdavis.health.comfort>

